	4010	3010	2010	1010	EasyAire
Weight Limit	250kg	200kg	200kg	180kg	180kg
Auto weight – fully weighs patient and adjusts pressure accordingly	✓	✓	*	×	×
Intelligent pressure sensing	✓	✓	✓	✓	*
Lower pressure and power failure alarms	✓	✓	✓	✓	✓
Panel locks	✓	✓	✓	×	✓
Switch on and everything else is automatically set up	✓	✓	✓	×	×
Static and dynamic therapy modes	✓	✓	✓	✓	*
Clinician override if required	✓	✓	✓	×	×
Autofirm - Gives static maximum inflation to allow nursing and care procedures and auto clinics off after 30 minutes	~	✓	×	×	×
Upright mode	✓	✓	×	×	×
Auto upright mode - Automatically optimises pressure when patient is sat up	✓	×	×	×	×
Autobody mass – adjusts pressure according to how the patients weight is distributed. i.e. A tall thin patient at 80kg would be given different pressure to a short wide patient also at 80kg	*	×	×	×	×
Auto switch off - switches off when the optimum internal pressure is reached	✓	✓	✓	✓	✓